Caution: Federal (U.S.A.) law restricts this device to sale by or on the order of a physician.

CooledTM Surgical Probe

Catalog Number: 600-004

U.S. Patents No. 5,651,780; 5,769,847; 5,797,905; 5,810,802; 5,871,523; 6,030,382; 6,106,522; 6,129,724; 6,387,092; 6,447,506; 6,471,699; 6,500,172; 6,142,994; 6,425,895; 6,120,496; 6,241,754; 6,416,505; 7,115,122; 6,830,576

Instructions for Use

TABLE OF CONTENTS

English	2
Dutch	6
French.....	10
German	14
Italian	18
Spanish	22

0086

Sterile – Single Use Only

ESTECH, Inc.
2603 Camino Ramon
Suite 100
San Ramon, CA 94583
USA
(925) 866-7111 Phone
(925) 866-7117 Fax
info@estech.com

European Authorized Representative:
Medical Partners BV
Groenendijk 39, Lage Zwaluwe
4926 RE Netherlands
+31.168.48 37 32 Phone
+31.168.48 47 90 Fax

ENGLISH

Instructions for Use

Carefully read all instructions prior to use. Observe all contraindications, warnings and precautions noted in these directions. Failure to do so may result in patient complications.

ESTECH relies on the physician to determine, assess and communicate to each patient all foreseeable risks of the procedure.

Description:

The ESTECH COBRA Cooled™ Surgical Probe (Probe) has a malleable distal section. These Probes are designed to conform to the specific anatomy of the tissue area to be coagulated.

The distal section of the probe has up to five 12.5 mm electrodes spaced 2 mm apart along the body. Any combination of coagulating electrodes may be used. The probe connects to the Cobra® Electrosurgical Unit (ESU) with a dedicated connecting cable. Directions for the use of the ESU and cable may be found in the ESU Operator's Manual.

For your convenience, a positioning tool is included in the packaging. This device is intended only to be used to help position and hold the probe against tissue.

Indications:

The Estech COBRA Cooled™ Surgical Probe is intended for the coagulation of cardiac tissue using radiofrequency (RF) energy during cardiac surgery. The System can be used during general surgery to coagulate soft tissues. The System may also be used to coagulate blood and soft tissue to produce hemostasis.

Contraindications:

- Local or systemic infection

Warnings:

- Contents supplied STERILE using an ethylene oxide (EO) process. Do not use if sterile barrier is damaged. If damage is found call your ESTECH representative.
- For single use only. Do not reuse, reprocess or resterilize. Reuse, reprocessing or resterilization may compromise the structural integrity of the device and/or lead to device failure which in turn may result in patient injury, illness or death. Reuse, reprocessing or resterilization may also create a risk of contamination of the device and/or cause patient infection or cross-infection, including, but not limited to, the transmission of infectious disease(s) from one patient to another. Contamination of the device may lead to injury, illness or death of the patient.
- The device(s) should be used by physicians thoroughly trained in the techniques of invasive surgical procedures.
- The operator should keep the set temperature and power limit as low as possible to achieve the desired end effect. This avoids excessive thermal damage to tissue, or collateral damage to adjacent tissue not intended for coagulation.
- Care should be taken to assure that the probe is not in contact with tissue other than that to be coagulated to avoid inadvertent tissue damage.
- Care should be taken when using the probe in proximity to vascular and nerve tissue to avoid inadvertent tissue damage.
- Care should be taken to thermally isolate the tissue to be coagulated when anatomically possible to avoid damage to unintended tissues or structures.
- Following RF coagulation, visual inspection of underlying tissues should be routinely performed to rule out the presence of inadvertent tissue damage.
- Care should be taken to assure that the probe is not in contact with other surgical instruments, staples or other objects while coagulating. Inadvertent contact with objects while coagulating could lead to conduction of RF Energy or heat and unintentional coagulation of tissues in contact with that object.
- Care should be taken when positioning surgical probe or positioning tool to prevent perforation or other damage to adjacent tissue.

**Precautions:**

- Before using, inspect for physical damage, including electrical insulation on the cables and the probe shaft. Replace damaged equipment.
- After use, dispose of product and packaging in accordance with hospital, administrative and/or local government policy.
- Do not attempt to operate the system before thoroughly reading the ESU Operator's Manual and Directions for Use. For future reference, keep these documents in a convenient, easily accessible place.
- While the distal portion of the probe is designed to be malleable and/ or flexible, to conform to the anatomy of the area to be coagulated, excessive or rough shaping of the probe may damage internal components of the probe.
- Pacemakers and implantable cardioverter/defibrillators can be adversely affected by RF signals. Refer to the manufacturer's Directions For Use.
- Dispersive Indifferent Patch (DIP) electrodes used with the ESTECH system should be applied carefully according to the manufacturer's directions. Poor or incomplete contact of the DIP electrodes may result in skin burns. The use of DIP electrodes which meet or exceed ANSI/AAMI 60601-2-2 requirements is recommended.
- The risk of igniting flammable gases or other materials is inherent in the application of RF power. Precautions must be taken to restrict flammable materials from the area where tissue coagulation is performed.
- Electromagnetic interference (EMI) produced by the ESU during the delivery of RF power may adversely affect the performance of other equipment.
- The Probe is intended for use only with the ESU, Instrument Cables and accessories.

Storage:

- Store in a cool dry place.

Inspection Prior to Use:

Carefully inspect the package prior to use for any breach of the sterile barrier or damage to the contents. If the sterile barrier integrity is compromised or the contents damaged, contact your ESTECH representative.

Directions for Use

1. Inspect the probe, positioning tool and all packaging materials carefully. Open the package using sterile technique.
2. Ensure that DIP electrodes are in good contact, according to the manufacturer's directions.
3. Connect the proximal end of the probe to the Instrument cable by plugging the cable into the port on the handle of the probe.
4. Gently shape the distal portion of the catheter to conform to the anatomy of the area to be coagulated. The positioning tool may be used to help position and hold the probe against the tissue.
5. Attach one female luer fitting to a fluid pump set at 500ml/hr flow rate. Either sterile water or normal saline can be used as a coolant. The other female luer fitting should be attached to a receptacle to catch coolant as it exits the probe. The fluid pump should be operating whenever the ESU is supplying energy.
6. Select electrodes to be activated on the ESU.

Ablation Time (seconds) COBRA Cooled™ Surgical Probe at 40W/Electrode*

Tissue Thickness	70° C	80° C
5mm	60	60
10mm	120	90

***Recommended ESU Settings: 40 Watts/coil, 70°C, 60 seconds (endocardial);
40 Watts/coil, 70°C, 120 seconds (epicardial)**

7. Activate selected electrodes by depressing the switch on the ESU.
8. Radiofrequency energy may be discontinued by depressing the switch on the ESU.

Complications

The following potential risks or discomforts may be associated with electrosurgical procedures. The frequency and severity of these events can vary, and may necessitate additional medical intervention, including surgery. Strict adherence to the forgoing instructions before use will help reduce the incidence of complications.

Allergic reaction, Arrhythmias, Cardiac or respiratory arrest, Cardiac valve damage, Chest pain, Damage to vessel intima or cardiac ultrastructures, Death, Hematoma / ecchymosis, Hemorrhage, Infarction, Infection, Perforation, Pericardial effusion, Pericarditis / pleuritis, Pseudoaneurysm, Pulmonary edema, Sinus or AV node injury, Stroke, Tamponade, Thrombosis, Vasovagal reaction.

Warranty and Limitations

ESTECH warrants that reasonable care has been used in the design and manufacture of this instrument. This warranty is in lieu of and excludes all other warranties not expressly set forth herein, whether express or implied by operation of law or otherwise, including, but not limited to, any implied warranties of merchantability or fitness for a particular purpose. Handling, storage, cleaning and sterilization of this instrument as well as other factors relating to the patient, diagnosis, treatment, surgical procedures, and other matters beyond ESTECH's control directly affect the instrument and the result obtained from its use. ESTECH's obligation under this warranty is limited to the repair or replacement of this instrument and ESTECH shall not be liable for any incidental or consequential loss, damage, or expense directly or indirectly arising from the use of this instrument. ESTECH neither assumes, nor authorizes any other person to assume for it, any other or additional liability or responsibility in connection with this instrument. ESTECH assumes no liability with respect to instruments reused, reprocessed or resterilized and makes no warranties, express or implied, including but not limited to merchantability or fitness for a particular purpose, with respect to such instrument.

For technical information contact:

ESTECH, Inc.

2603 Camino Ramon, Suite 100

San Ramon, CA 94583 USA

Phone: 925-866-7111

Fax: 925-866-7117

email: info@estech.com

European Authorized Representative

Medical Partners BV

Groenendijk 39

Lage Zwaluwe

4926 RE Netherlands

Phone: +31.168.48 37 32

Fax: +31.168.48 47 90

THIS PRODUCT CONTAINS NO LATEX.

THIS PRODUCT MAY CONTAIN DEHP.

Graphic Symbols for this Device Labeling

<p>Catalogue Number</p> 	<p>Caution-See instructions for Use</p> 	<p>Manufacturer/Company Address</p> 	<p>Sterilized using Ethylene Oxide</p> 
<p>Batch Code</p> 	<p>Single Use Only</p> 	<p>Expiration Date</p> 	<p>Authorized Representative in the European Union</p> 
<p>Quantity</p> 	<p>Date of Manufacture</p> 	<p>Non-Pyrogenic</p> 	<p>Does not contain Natural Rubber Latex</p> 
<p>Contains no di (2-ethylhexyl) phthalate (DEHP) released from polyvinyl chloride (PVC)</p> 	<p>"Consult Instructions For Use"</p> 	<p>CE Mark with Notified Body - class II and III medical devices comply with EC Directive 93/42/EEC</p> 	<p>"Caution: U.S. Federal law restricts this device to sale by or on the order of a physician or other licensed practitioner."</p> <p>Rx ONLY</p>

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